



A free toolkit for UK clinicians
Everything but the ambulance in your pocket

- ✓ Tap to Count Respiratory Rate
- ✓ Observations & Handovers
- ✓ NEWS2, GCS & Frailty Scoring
- ✓ Cardiac Arrest Logging
- ✓ Mental State Exam prompt
- ✓ Timer system

Plus many more tools and clinical aides

11:51

SnapMedic

← Cardiac Arrest

Mode: ALS Switch to BLS

110 BPM

Stop Metronome

Actions Log

Arrest Length 0m 26s Adrenaline 0m 22s Rhythm Check 0m 20s

Interventions

Shock (VF) Shock (VT)

Rhythm Check (PEA) Rhythm Check (Asystole)

2nd Adrenaline 1st Amiodarone

Advanced Airway ROSC

Custom note... + Add

12:11

SnapMedic

← Observations

VITAL SIGNS 3m ago

Blood Pressure 165 / 94 mmHg Heart Rate 114 bpm

Respiratory Rate 35 /min SpO₂ 97 %

ETCO₂ 7.2 kPa Blood Sugar 5.5 mmol/L

Temperature 39.5 °C Pain Score 0 /10

GLASGOW COMA SCALE 3m ago

Eye Verbal Motor Total

4 5 6 15

PUPIL ASSESSMENT 3m ago

Left Pupil Right Pupil

Size (mm) Size (mm)

3 mm 3 mm

Reactivity

12:00

SnapMedic

← Respiration Counter

Breaths 6 Duration 24.6s

Breathing Pattern

Breathing intervals (taller = longer pause between breaths)

Regular Variable Irregular

Pattern: Variable · Avg: 4.9s

Reset Counter

➦ Add to Observations

15 breaths/min



A free toolkit for UK clinicians

Everything but the ambulance in
your pocket



What are you waiting for?

snapmedic.app